## **AED CABINET ALARMS**







## STROBE ALARM (14SA)

Acts in conjunction with the standard Commander Alarm to provide a visual strobe light and a 110dB warning horn.

- Flashes strobe light when door is opened.
- 110dB warning horn sound alert.
- On cabinets with at least 3" or wider trim, the strobe alarm will be installed. For cabinets under 3" wide trim, the alarm will be shipped loose for installation on the wall.

Note: This alarm only operates when connected with the Commander Alarm.

## **COMMANDER ALARM (14A)**

## Standard on all AED cabinets

- · Installs on inside of cabinet.
- Plunger switch is activated when the cabinet door is opened.
- 85dB warning horn sound alert.
- Alarm is deactivated when the door closes or the key switch is turned to the off position.
- Keyed on/off switch for AED cabinets only.
- 9-Volt battery not included. (Part # 14-9V)

It is the responsibility of the customer to communicate specifications and product performance requirements that are needed for their application. This document is intended to convey potential product configurations and offerings. It is not to be used to imply or represent a specification or rating, of finished products shipped. Detailed, specific submittals can be created upon request. It is the responsibility of the customer to specify product accordingly to comply with local laws, codes, and ordinances.

